

ZSI 375 PF ARTIFICIAL URINARY SPHINCTER

First Artificial Urinary Sphincter Pre-Filled



- Pre-connected / Pre-filled / 2 components device.
- One size fits all. Adjustable cuff from 3.75cm to 5cm (circumference).
- Ergonomic cuff moulded in Curve forme.
- Dissection-Insertion. No preparation.
- Reduced risk of infection due to short time in the OR.
- Possibility of increasing issued pressure.


Pre filled



Possibility of increasing
Issued pressure



1 3.75 to 5 cm adjustable
Cuff, moulded in a curve



2 Deactivated when the button
(▼) portrudes at opposite
Side of the pump button



3 Pump button to empty
The cuff (press-release)



4 Activated when the button (▼)
portrudes at the side of the
Pump button



ZSI 475 HYDRAULIC PENILE IMPLANT




- All Standard sizes come in one box.
- Accessories for implantation are provided in the box / except dilators.
- Reduce time of kit preparation.
- PVP coating.
- All components are biocompatible.
- Larger diameter when inflated, resistive multiple folding, better rigidity.
- Extra large size 22cm available. *




ZSI 100 MALLEABLE PENILE IMPLANT

- Pair of implants adjustable through 12-25cm.
- Adjustable distally/proximally for easy implantation.
- Flexible distal section discreetly conceals device.
- Flexible distal section prevents gland discomfort in trousers.
- High quality memory form for optimal comfort both erect and resting.




ZSI 100 FTM MALLEABLE PENILE IMPLANT FOR PHALLOPLASTY




DISTAL PART
Adjustable by cutting 5mm by 5mm.

CENTRAL PART
Cylinder of Ø 22mm in width. Designed for comfort and pleasant appearance.

PROXIMAL PART
30mm by 25mm. Is fixed securely onto the pubic bone.



ZSI 475 FTM HYDRAULIC PENILE IMPLANT FOR PHALLOPLASTY



DISTAL PART
glands diameter 25 mm.

CENTRAL PART
main body diameter Ø 21-22 mm.

PROXIMAL PART
30mm by 25mm.

